



# Neater Steady Double Straw (ND-KV2)

The Neater Steady Double Straw provides a stable and adjustable holder for two flexible straws. A non-return valve holds the liquid in the straw to make sucking easier.




## To assemble the Neater Steady Double Straw:

- 1 Clamp to table.
- 2 Clip in flexible straw.
- 3 Bend to a convenient position. Be careful to AVOID sharp bends and DO NOT pull the bendy arm hard.
- 4 Use Velcro to tidy the set up .



The straw can be pre-filled using a syringe.



To help with the first suck - the straw can be pre-filled using a syringe.

With the plunger of the syringe pushed in, insert the syringe nozzle into the top end of the straw and gradually pull out the plunger to suck the drink up the straw.

You only need to suck gently to use the Neater Double Steady Straw

The non-return valve at the bottom of the rigid straw prevents the drink flowing back into the mug. If your suck is very gentle, then it will take several sucks to fill the straw. Once it is full though, every suck will give you some drink.

If your suck is very weak, you could try shortening the straw and placing the mug in a higher position (though never above the level of your mouth).

Some people may benefit from a short length of the rigid tubing (The Hard mouth piece as depicted on the front) being pushed all the way inside the mouthpiece end of the flexible straw.



Using the Straw Cleaning Brush (ND-SC)

To use, feed the brush handle end of the cleaning brush (which can be found in an envelope in your ND-KV2 pack) without the bristles all the way into the straw and pull through. You can repeat this as necessary.

Food safe cleaning and/or sterilising solution can be used. Rinse the straw once you are finished.

## Safety

- The Neater Drinker should only be used for drinking consumable liquids.
- It should not be used by people with swallowing difficulties or who are liable to choke, except under the supervision of a trained Speech and Language Therapist.
- Drinks should not be scalding hot.
- Your mouth should always be higher than the mug or siphoning can occur and cause choking.
- Clean all parts before use and check that no cracks have appeared in any components.
- The syringe is for filling the straw only.
- Take care not to stick the straw in your eye or into the back of your throat.
- Keep the stop piece on the straw to help ensure its visibility and to help you avoid putting the straw too far back in your mouth.
- Soft mouthpiece tubing may reduce injury.

For other drinker options, please visit our website!



**Neater.co.uk**  
Independence • Dignity • Choice

info@neater.co.uk

+ 441298 23883